


























Oticon Hearing Instrument Portfolio

Invisible In Canal	In Canal		In The Ear	Receiver In The Ear			Behind The Ear			
IIC Invisible In Canal	CIC Completely In The Canal	ITC In The Canal	ITE In The Ear	miniRITE Mini Receiver In The Ear	miniRITE T Receiver In The Ear	miniRITE R Receiver In The Ear	BTE Behind The Ear	Plus Power Power Behind The Ear	BTE 13SP SuperPower Behind The Ear	BTE 13UP Ultra Power Behind The Ear
										
			 Full Shell (pictured above)  Half Shell				  Available in both thin tube and hook solution. Only available for Ruby and Siya	  Available in both thin tube and hook solution.		
<p>Invisible In Canal</p> <p>Suitable for first-time and experienced users with a mild to severe hearing loss, who are looking for the ultimate discreet hearing solution. This style sits deepest within the ear canal so it is not visible at all.</p>	<p>Completely In The Canal</p> <p>Suitable for mild to severe hearing loss. This tiny custom style rests deep in the ear canal.</p>	<p>In The Canal</p> <p>Suitable for mild to profound hearing loss. This style rests in the ear canal.</p>	<p>In The Ear</p> <p>Suitable for mild to profound hearing loss. The FS style fills the concha of the ear while the HS style fills less than half of the concha.</p>	<p>Mini Receiver In The Ear</p> <p>Suitable for mild to profound hearing loss. A very small attractive RITE style that fits discreetly behind the ear transmitting sound via a thin wire.</p>	<p>Mini Receiver In The Ear</p> <p>Suitable for mild to profound hearing loss. MiniRITE T is a discreet style, based on the popular miniRITE, and features both telecoil and a convenient double pushbutton for easy volume and program control.</p>	<p>Mini Receiver In The Ear</p> <p>Suitable for mild to profound hearing loss. This is a discreet style powered by rechargeable lithium-ion battery. MiniRITE R features telecoil and a convenient double pushbutton.</p>	<p>Behind The Ear</p> <p>Suitable for mild to severe hearing loss. Available in both thin tube and hook solutions, this style rests behind the ear and transmits sound to the ear via a tube.</p>	<p>Plus Power BTE</p> <p>Suitable for severe to profound hearing loss. Available in both thin tube and hook solutions, this style rests behind the ear and offers more power than the BTE in a discreet housing.</p>	<p>BTE 13SP</p> <p>Suitable for severe to profound hearing loss. Super power hearing aid with a 13 battery. It has separate push buttons for programs and volume for easy usage and control. It features T-coil, optional LED indications and FM support.</p>	<p>BTE 13UP</p> <p>Suitable for severe to profound hearing loss. Ultra power hearing aid with a 675 battery. This style has separate push buttons for programs and volume for easy usage and control. It features T-coil, optional LED indications and FM support.</p>

705-360-6680

Timmins and District Hospital
700 Ross Avenue East, Timmins, ON P4N 8P2

ethierhearingclinic.com

